

3C\TEMP-USB -200C

- Multi-use
- Cryogenic temperatures (liquid nitrogen)
- Dual sensors
- Real-time LCD display
- User-selectable settings
- Easy to use. Just push the start button to begin recording.
- Download all recorded data to your PC with a compatible cable and Marathon's MaxiThermal2 software.



Internal Sensor: -29°C to +72°C
External Probe: -200°C to +72°C



幸託有限公司
XIN TOP CORPORATION

TEL : (02)2598-1199
FAX : (02)2596-2331

E-mail : info@xintop.com
Website : www.xintop.com